

soluCLEAN

Daily Toilet Cleaner




Dilution ratio

1 sachet per 1L bottle

This Daily toilet cleaner is for internal surfaces of sanitary ware, squeeze the bottle to apply liquid under the rim of the toilet bowl. Reduces odour leaving the toilet clean and fresh. Squirt around the toilet, leave to stand for 5 minutes, scrub with a brush and flush away.



Instructions




1. Place the sachet in the bottle.



2. Half fill with warm water.



3. Screw the bottle head back on and shake the bottle vigorously until sachet fully dissolves. Leave to stand for 15 minutes.



4. Top up with warm water to fill line. Shake and ready to use.

Solupak Product Code:

SCSP1000TCP-10

Solupak Bottle code:

10 pack- SCPPBOT1000TCP

Best practice

We don't like harming people or the environment so our products are not hazardous when they are diluted. The powder in the sachet is hazardous however we have wrapped it up safely in water soluble film so you don't come into contact with it.

Please make sure you:

Store the sachets in a dry area and keep them inside the packet until you need them.

Only handle the sachets with dry hands. If any sachets are damaged please do not handle them without wearing PPE and always avoid the powder making contact with your eyes or skin.

Read and follow the safety data sheet provided.

Always follow the colour coding.

Once diluted the solution is not classified with any hazard.

Where to use



Urinal



Toilet



Housekeeping

Classification

POWDER ONLY

ONCE DILUTED



NOT CLASSIFIED UNDER CLP


But the sachets outer film keeps you safe

Type:

Squirt

pH in use:

3-4



Colour:

Red