

# soluCLEAN

Sustainably Derived

## Disinfectant


### Dilution ratio

1 sachet per 750ml trigger bottle

Independently tested to EN 13697, also conforms to BSEN 1276, PAS 2424 Bacteriostatic. (list of materials above). Kills odours and Leaves behind long lasting mango fragrance. No pre-clean required- cleans and disinfects in one use.



## Instructions



1. Add water to the fill line.
2. Place the sachet in the bottle.
3. Screw the spray head back on and shake the bottle.
4. Once the sachet has dissolved, spray onto the surface, leave for 5 minutes, then wipe.

## Where to use



Germ Kill



Work Surface

Solupak Product Code: | DCSP750DDP-10

Solupak Bottle Code: | SCPPBOTTLEDDP

## Classification

POWDER ONLY

ONCE DILUTED



NOT CLASSIFIED UNDER CLP

## Best practice

We don't like harming people or the environment so our products are not hazardous when they are diluted. The powder in the sachet is hazardous however we have wrapped it up safely in water soluble film so you don't come into contact with it.

Please make sure you:

Store the sachets in a dry area and keep them inside the packet until you need them.

Only handle the sachets with dry hands. If any sachets are damaged please do not handle them without wearing PPE and always avoid the powder making contact with your eyes or skin.

Read and follow the safety data sheet provided.

Always follow the colour coding.

Once diluted the solution is not classified as with any hazard.

Type: | TRIGGER

pH in use: | 6.5



Acidic | Neutral | Alkaline

Colour: | Blue