soluCLEAN Sustainably Derived

Washroom Disinfectant



I sachet per 750ml bottle

Indpendently tested to EN 13697, also conforms to BSEN 1276, PAS 2424 Bacteriostatic. (list of materials above). Kills odours and Leaves behind long lasting mango fragrance. No pre-clean required- cleans and disinfects in one use.



Instructions Where to use GO, once the sachet has SHAKE, screw FILL to ADD the sachet the spray head back on and dissolved, spray onto the the wate surface, leave for 5 minutes and then wipe shake the bottle Classification Solupak SCFS10WDP-10 **Product Code:** ONCE DILUTED POWDER ONLY NOT CLASSIFIED UNDER CLP Solupak SCPPBOTTLE750WDP **Bottle Code:** But the sachets outer film keeps you safe **Best practice** TRIGGER Type: We don't like harming people or the environment so our products are not hazardous when they are diluted. The powder in the sachet is hazardous however we have wrapped it up safely in water soluble film so you don't come into contact with it. pH in use: 2.3 Please make sure vou: Store the sachets in a dry area and keep them inside the packet until you need them. Only handle the sachets with dry hands. If any sachets are damaged please do not 2 3 4 5 6 7 8 9 10 11 12 13 14 0 1 handle them without wearing PPE and always avoid the powder making contact with your eyes or skin. 4 Read and follow the safety data sheet provided. Acidic Neutral Alkaline Always follow the colour coding. Once diluted the solution is not calssified as with any hazard. **Colour:** Red